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

2008311 // BIG BLUE WAVE

nitrile disposable gloves

PPE risk category III | Medical device class I

EN

Manufacturer's instructions and information

Information brochure for personal protective equipment (PPE) according to Regulation (EU) 2016/425, annex II point 1.4. Please read this information brochure carefully before using the PPE. You are obligated to enclose this information brochure when passing on the PPE or to hand it over to the recipient of the PPE. For this purpose, this information brochure may be reproduced without restriction.

nitrile disposable gloves	PPE risk category III Medical device class I
Size(s):	S-XL
Colour code:	2200
Certification:	EN ISO 21420, EN ISO 374, EN 421

The CE marking certifies that the product complies with the essential health and safety requirements of Regulation (EU) 2016/425. The EU Declaration of Conformity can be viewed at [doc.nitras-medical.com](#).

This product is personal protective equipment (PPE) classified as Risk Category III. It protects you against risks that could lead to very serious consequences, such as death or irreversible damage to health. This product offers protection against: chemicals, microorganisms, viruses, and radioactive contamination from solid particles. Any uses other than those specified above are expressly excluded. This product therefore does not, amongst other things, provide protection against: mechanical risks, cold, thermal risks (heat and/or fire), electric shock, radiation, or work involving high-pressure jets. Please observe the pictograms, instructions and corresponding performance levels provided.

Storage / Use / Inspection: Store in a cool, dry place. Keep away from direct sunlight, UV rays and ozone sources. Do not store in a bent position or under heavy loads. Store and transport the product in its original packaging if possible. Factors such as light, moisture, temperature and natural material changes over a longer period of time can result in a change in the product properties. It is not possible to provide exact information on the storage time and service life of PPE, as both parameters depend, among other things, on the type of storage, temperature, humidity, degree of wear and intensity of use. Therefore, check this product for damage or material changes (e.g. brittle, cracked coatings/materials, holes, colour changes, etc.) after prolonged storage and before and after each use. Check this product before each use to ensure it is suitable for the intended activity and is the correct size. Unsuitable or defective products must be disposed of and must not be used under any circumstances. The size of the product may differ from the specifications, e.g. due to stretching.

All performance data was determined through testing under laboratory conditions. It is therefore recommended to check whether the PPE is suitable for the intended use, as conditions in the workplace may differ from those of the type examination depending on various parameters (e.g. temperature, abrasion, intensity of use). If PPE has already been used, it may offer lower performance due to the degree of wear. The manufacturer accepts no responsibility for improper use of the product.

Instructions for wearing the item: Ensure that your hands are clean and dry before putting on gloves. Insert your fingers into the respective gloves and pull the gloves loosely over your hands by the knitted cuffs or cuffs. Make sure they fit correctly. Gloves should fit snugly and tightly on the palm, fingers and between the fingers. Fingernails, jewellery and excessive stretching and pulling can damage the gloves. After use, gloves should be removed in such a way that the outside does not come into contact with clothing or skin, as these may be visibly or invisibly contaminated with harmful substances. Gloves should therefore be removed with the inside facing outwards. To do this, first loosen the fingertips of the glove from your fingers. The knitted cuff or cuff can then be turned inside out to remove the glove. To maintain the comfort of the glove, it should be cleaned after each use in accordance with the cleaning and maintenance instructions. Depending on requirements, this can and should be done while wearing the gloves.

Before starting work (after breaks and, if necessary, after washing hands), a suitable skin protection product can be used. During work (before breaks and before finishing work), a suitable skin cleanser can be used. After work (after the last hand wash), a suitable skin care product can be used.

Cleaning / Maintenance: This product is intended for single use only and must be disposed of regularly. Expiry date: see packaging.

Disposal: Dispose of this product with household waste. After intentional or unintentional contact with chemicals, this product may be contaminated with environmentally harmful or hazardous substances. In this case, disposal must be carried out in accordance with local regulations.

Special notes: PPE may cause allergic reactions in sensitive individuals. Special caution is advised in cases of known hypersensitivity.

General explanations of performance levels achieved: 1-6 / A-F: test result achieved (the higher, the better) 0: minimum performance level not achieved X: not tested or not applicable due to the material or design All tests were carried out under laboratory conditions on the palm of the hand, and the respective performance levels were determined on this basis.

Classes and breakthrough time (minutes)

Class 1: > 10 minutes / Class 4: > 120 minutes

Class 2: > 30 minutes / Class 5: > 240 minutes

Class 3: > 60 minutes / Class 6: > 480 minutes

Results in accordance with EN ISO 374-2:2019.7.2 / 7.3: passed

This information does not provide details regarding the actual duration of protection in the workplace or the distinction between mixtures and pure chemicals. Chemical resistance was assessed under laboratory conditions on samples taken solely from the palm of the hand (except where the glove is 400 mm or longer – in which case the cuff is also tested) and relates exclusively to the chemicals tested. It may differ if the chemical is used in a mixture. It is recommended that you check whether the gloves are suitable for the intended use, as workplace conditions may differ from those of the type test depending on temperature, abrasion and degradation. If protective gloves have already been used, they may offer reduced resistance to hazardous chemicals due to changes in their physical properties. Degradation caused by contact with chemicals, movement, snagging, friction, etc., can significantly reduce the actual service life. In the case of aggressive chemicals, degradation may be the most important factor to consider when selecting chemical-resistant gloves. Before use, the gloves must be checked for any faults or defects. Decontamination of chemical and biological contamination must be carried out in a specific manner. The nature of the contamination must be known both qualitatively and quantitatively in order to assess the degree of decontamination. In any type of decontamination, personal protection is essential to prevent harm to the individual and the environment. This means that, along with the contaminants, the decontamination agents and personal protective equipment (water, cleaning agents, brushes, filters, gloves and clothing) must be collected and either disposed of properly or cleaned in a specific manner. In principle, personal protective equipment should be removed and set aside in such a way that the outer surface does not come into contact with clothing or skin. Protective gloves should therefore be removed so that the inside faces outwards. The decontamination instructions apply exclusively to gloves that can be reused. Other gloves are not suitable for decontamination and are intended for single use only. Please refer to the information under Cleaning / Maintenance.

Test parameter	Degradation (%)
J - n-Heptane	59
K - Sodium hydroxide 40%	-29,4
T - Formaldehyde 37%	22,7

EN ISO 374-5:2016 Protective gloves against dangerous chemicals and micro-organisms			
ISO 374-5:2016	Test parameter	Performance levels	Test result
 VIRUS	Penetration	-	-
	Protection against viruses	-	-

These gloves protect against microorganisms (bacteria and fungi). If the label 'VIRUS' appears below the pictogram, the gloves also protect against viruses. If this label is not present, the gloves have not been tested for protection against viruses. Resistance to penetration has been assessed under laboratory conditions and applies exclusively to the samples tested.

EN ISO 21420:2020 + A1:2024 Protective gloves – General requirements and test methods			
EN ISO 21420	Test parameter	Performance levels	Test result

If there is a risk of getting caught in moving machine parts, gloves must not be worn.

EN 421:2010 Protective gloves against ionizing radiation and radioactive contamination			
EN 421:2010	Test parameter	Performance levels	Test result
			

Manufacturer

Year and month of manufacture

Recycling symbol (only for FR)

Read the manufacturer's instructions and information

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