# 8305 // PROTECT 300

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EN

#### Manufacturer's instructions and information

Information brochure for personal protective equipment (PPE) according to Regulation (EU) 2016/425, annex II point 1.4. Please read this information brochure carefully before using the PPE. You are obligated to enclose this information brochure when passing on the PPE or to hand it over to the recipient of the PPE. For this purpose, this information brochure may be reproduced without restriction.

Protective gloves	Risk category III
Size(s)	S-XL
Certification	EN ISO 374
Notified body	ANCCP Certification Agency Srl Via Dello Struggino, 6 57121 Livorno Italy 0302

The CE marking certifies that the product complies with the essential health and safety requirements of Regulation (EU) 2016/425. The EU declaration of conformity can be viewed at [www.doc.nitras.de](http://www.doc.nitras.de). This product is personal protective equipment of risk category III. It protects you against risks that can lead to very serious consequences such as death or irreversible damage to health. This product offers protection against: chemicals, microorganisms. Other areas of application than those mentioned above are expressly excluded. This product therefore provides, among other things, no protection against: mechanical hazards, cold, thermal risks (heat and/or fire), electric shock, radiation, high-pressure jets. Please note the pictograms, notes and the corresponding performance levels.

Storage / use / servicing: Store in a cool, dry place. Keep away from direct sunlight, UV rays or ozone sources. Do not store in buckled condition or under weight load. If possible, store or transport the product in its original packaging. Influences such as light, humidity, temperature and natural changes in materials over a longer period of time can lead to changes in product properties. Exact information on storage time and service life of the PPE is not possible, since both parameters depend on the respective type of storage, temperature, humidity, degree of wear and intensity of use, among other things. Check this product for damage or material changes (e.g. brittle, cracked coatings / materials, holes, colour changes etc.) after prolonged storage and before and after each use. Before each use, check this product for suitability for the intended activity and for the correct size. Unsuitable or defective products must be disposed of and never used. The size of the product may differ from the specifications, e.g. due to stretching.

All performances were determined by tests under laboratory conditions. It is therefore recommended to check whether the PPE is suitable for the intended use, as the conditions at the workplace can differ from those of the type examination depending on various parameters (e.g. temperature, abrasion, intensity of use). If PPE has already been used, it can offer lower performance due to the degree of wear. The manufacturer accepts no responsibility for any improper use of the product. Instructions for wearing the product: Make sure your hands are clean and dry before putting on gloves. Insert your fingers into the respective glove and pull the glove loosely over your hand on the knitted wrist or cuff. Make sure that the fit is correct. Gloves should have a tight and snug fit on the palm of the hand, fingers and gaps between fingers. Fingermats, jewellery, excessive stretching and pulling can damage the gloves. Gloves should be taken off after use in such a way that the outside of the gloves does not come into contact with clothing or skin, as the gloves can be visibly and invisibly contaminated with harmful substances. Accordingly the inside must come outwards. First remove the fingertips of the glove from your fingers. The knitted wrist or cuff can then be rolled outwards in order to remove the glove. To ensure that the glove retains its comfort, it should be cleaned after each use in accordance with the cleaning and maintenance instructions. If necessary, this can and should be done while wearing the gloves. A suitable skin protection product can be used before starting work (after breaks and if necessary after washing the hands). During work (before breaks and before end of work) a suitable skin cleanser can be used. After work (after the last washing of the hands) a suitable skin care product can be used. Cleaning / maintenance: This product is intended for single use and must be disposed of regularly. Expiry date: See packaging.

Disposal: Dispose of with household waste. This product may be contaminated by environmentally harmful or hazardous substances after intended or unintended contact with chemicals. In this case, disposal must be carried out in accordance with the local legal regulations. Special notes: PPE can cause allergic reactions. Special care is recommended in case of known hypersensitivity.

General explanations of achieved performance levels  
1-6 Achieved test result (the higher, the better)  
0 Minimum performance level not achieved  
X Not tested or not applicable due to the material or design  
All tests were carried out under laboratory conditions on the palm of the hand. Respective performance levels were determined on this basis.

Test parameter	Performance level	Test result
Dexterity	1-5	5

This article is intended to protect the user's hand from substances, liquids or chemicals for a short period of time. Therefore, these are gloves for special applications and the length of the gloves deviates from the dimensions according to EN 420. If protection is also required on the wrist, other gloves have to be used. If there is a risk of getting caught in moving machine parts, gloves must not be worn.

#### EN 420:2003 + A1:2009 Protective gloves - General requirements and test methods

Test parameter	Performance level	Test result
Dexterity	1-5	5

#### EN ISO 374-1:2016 Protective gloves against dangerous chemicals and micro-organisms

ISO 374-1:2016/Type C	Test chemical	Code letter	Class	Test result
	n-heptane	J	1-6	1
	Sodium hydroxide 40%	K	1-6	6
	Formaldehyde 37%	T	1-6	4

Class	Breakthrough time (minutes)	Class	Breakthrough time (minutes)
1	> 10	4	> 120
2	> 30	5	> 240
3	> 60	6	> 480

Results according to EN 374-2:2014, 7.2 / 7.3: Passed  
Results according to EN 374-4:2013:

Test chemical	Degradation (%)
n-heptane	60,8
Sodium hydroxide 40%	-29,4
Formaldehyde 37%	22,7

#### EN ISO 374-5:2016 Protective gloves against dangerous chemicals and micro-organisms

ISO 374-5:2016
VIRUS

This information does not indicate the actual duration of protection at the workplace and the distinction between mixtures and pure chemicals. The resistance to chemicals has been assessed under laboratory conditions on samples taken only from the palm of the hand (except where the glove is 400 mm or longer - in which case the cuff is also tested) and refers exclusively to the chemicals tested. It can be different if the chemical is used in a mixture. It is recommended to check whether the gloves are suitable for the intended use, as the working conditions at the workplace may differ from those of the type test depending on temperature, abrasion and degradation. If protective gloves have already been used, they may be less resistant to hazardous chemicals due to changes in their physical properties. Degradation, movement, thread pulling, friction etc. caused by contact with chemicals can considerably reduce the actual application time. For aggressive chemicals, degradation may be the most important factor to consider when selecting chemical resistant gloves. Before use, the gloves must be checked for any faults or defects.

The decontamination of chemical and biological contamination must be carried out specifically. The contamination must be known both qualitatively and quantitatively in order to be able to make a statement about the degree of decontamination. In any type of decontamination, self-protection is important to avoid an endangerment of the person and the environment. This means that, together with the contaminants, the materials used for decontamination and personal protective equipment (water, cleaning agents, brushes, filters, gloves and clothing) must be collected, disposed of or specifically cleaned. In principle, personal protective equipment should be taken off and put away in such a way that the outside of the protective equipment does not come into contact with clothing or skin. Protective gloves must therefore be removed in such a way that the inside of the glove comes outwards.

These gloves protect against microorganisms (bacteria and fungi). The resistance to penetration was assessed under laboratory conditions and relates exclusively to the samples tested. Not checked against viruses.

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