

TECHNICAL DATASHEET

004332 // SOFT PROTECT PLUS

Material: Non-woven fabric, dark brown
Others: One size

Packaging Unit: 500 pieces



Description

NITRAS SOFT PROTECT PLUS, medical face mask, made of very soft and fibreglass-free non-woven fabric, 4-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 Type IIR, bacterial filtration efficiency > 98%, 10 boxes, 50 pieces each

Construction

Layer 1: coloured polypropylene layer
Layer 2: coloured polypropylene layer
Layer 3: white melt blown filter
Layer 4: white polypropylene layer

Two elastic, latex-free elastic bands

Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

Technical specification

 $\begin{array}{lll} \mbox{Layer 1:} & 23 \ g/sqm \pm 2 \ g/sqm \\ \mbox{Layer 2:} & 23 \ g/sqm \pm 2 \ g/sqm \\ \mbox{Layer 3:} & 20 \ g/sqm \pm 2 \ g/sqm \\ \mbox{Layer 4:} & 20 \ g/sqm \pm 2 \ g/sqm \\ \mbox{Measurement:} & 17.5 \ cm \ x \ 9.5 \ cm \pm 2 \ mm \\ \end{array}$

Weight: $5 g \pm 0.2 g$

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural materical changes during a period could implicate change in the properties.