

004335 // SOFT PROTECT PLUS

Face masks
Medical device class I



Breathable whilst providing reliable protection
Integrated nosepiece for optimal fit
Latex-free ear loops for optimal hold

004335 // SOFT PROTECT PLUS

Material: Polypropylene, fossil

Miscellaneous: 4-ply

Packaging unit: 10 piece(s)

NITRAS SOFT PROTECT PLUS, medical face mask, taupe (colour code: 5500), non-woven fabric, 4-ply, fiberglass free, breathable whilst providing reliable protection, very soft, integrated nosepiece for optimal fit, round and latex-free ear loops for optimal hold, colour-coordinated, EN 14683 Type IIR, bacterial filtration efficiency > 98%, box per 50 pieces

EN 14683
Type IIR

Industries

-  Beauty Edition
-  Care Edition
-  Dental Edition
-  Medical Edition

Advantages

- + Breathable whilst providing reliable protection
- + Integrated nosepiece for optimal fit
- + Latex-free ear loops for optimal hold



Further images



Colour variants



004335 // SOFT PROTECT PLUS
Colour code: 5500



004331 // SOFT PROTECT PLUS
Colour code: 2100



004333 // SOFT PROTECT PLUS
Colour code: 6500



004334 // SOFT PROTECT PLUS
Colour code: 6300



004336 // SOFT PROTECT PLUS
Colour code: 6600



004338 // SOFT PROTECT PLUS
Colour code: 1000



004339 // SOFT PROTECT PLUS
Colour code: 6400

Additional information

	Material	Colour	Grammage
Shell fabric	100% Polypropylene	taupe	

Physical properties

Grammage	23 g/qm
Weight	5 g
Dimensions	17,5 cm x 9,5 cm ± 3 mm

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904166240	4059904166257	0,52 m	0,2 m	0,26 m	3,09 kg	56

Customs tariff number: 63079095110

Country of origin: China