

TECHNICAL DATASHEET

004336 // SOFT PROTECT PLUS

Material:	Polypropylene, bordeaux
Others:	One size
Packaging Unit:	500 pieces



Description

NITRAS SOFT PROTECT PLUS, medical face mask, made of very soft and fibreglass-free non-woven fabric, 4-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 Type IIR, bacterial filtration efficiency > 98%, 10 boxes, 50 pieces each

Construction

Layer 1: coloured polypropylene layer
 Layer 2: coloured polypropylene layer
 Layer 3: white melt blown filter
 Layer 4: white polypropylene layer
 Two elastic, latex-free elastic bands
 Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

Technical specification

Layer 1: 23 g/sqm ± 2 g/sqm
 Layer 2: 23 g/sqm ± 2 g/sqm
 Layer 3: 20 g/sqm ± 2 g/sqm
 Layer 4: 20 g/sqm ± 2 g/sqm
 Measurement: 17.5 cm x 9.5 cm ± 2 mm
 Weight: 5 g ± 0.2 g

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural material changes during a period could implicate change in the properties.