

# TECHNICAL DATASHEET

#### 004338 // SOFT PROTECT PLUS

Material: Polypropylene, black
Others: One size
Packaging Unit: 500 pieces



#### **Description**

NITRAS SOFT PROTECT PLUS, medical face mask, made of very soft and fibreglass-free non-woven fabric, 4-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 Type IIR, bacterial filtration efficiency > 98%, 10 boxes, 50 pieces each

### Construction

Layer 1: coloured polypropylene layer
Layer 2: coloured polypropylene layer
Layer 3: white melt blown filter
Layer 4: white polypropylene layer

Two elastic, latex-free elastic bands

Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

## **Technical specification**

Layer 1:  $28 \text{ g/sqm} \pm 2 \text{ g/sqm}$ Layer 2:  $28 \text{ g/sqm} \pm 2 \text{ g/sqm}$ Layer 3:  $20 \text{ g/sqm} \pm 2 \text{ g/sqm}$ Layer 4:  $20 \text{ g/sqm} \pm 2 \text{ g/sqm}$ Measurement:  $17.5 \text{ cm } \times 9.5 \text{ cm} \pm 2 \text{ mm}$ 

Weight:  $5 g \pm 0.2 g$ 

#### Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural materical changes during a period could implicate change in the properties.