

TECHNICAL DATASHEET

004339 // SOFT PROTECT PLUS

Material: Polypropylene, pink
Others: One size
Packaging Unit: 500 pieces



Description

NITRAS SOFT PROTECT PLUS, medical face mask, made of very soft and fibreglass-free non-woven fabric, 4-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 Type IIR, bacterial filtration efficiency > 98%, 10 boxes, 50 pieces each

Construction

Layer 1: coloured polypropylene layer
Layer 2: white melt blown filter
Layer 3: white polypropylene layer
Layer 4: white polypropylene layer

Two elastic, latex-free elastic bands

Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

Technical specification

Layer 1: 23 g/sqm \pm 2 g/sqm Layer 2: 25 g/sqm \pm 2 g/sqm Layer 3: 23 g/sqm \pm 2 g/sqm Layer 4: 40 g/sqm \pm 2 g/sqm Measurement: 17.5 cm x 9.5 cm \pm 2 mm

Weight: $5 g \pm 0.2 g$

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural materical changes during a period could implicate change in the properties.