

004546

PPE risk category III



Breathable and water-repellent
Perfect hold
Protection against infection agents

004546

Material: Polypropylene, light blue

Miscellaneous: Size 2XL

Packaging unit: 50 piece(s)

NITRAS medical gown, light blue (colour code: 2300), polyethylene with laminate film, durable, breathable and water-repellent, length approx. 110 cm, pleasant wearing comfort, with round collar, soft cotton ribbed cuffs for a better fit and perfect hold, to tie behind the neck and body, protection against infection agents, Type PB [6]-B, single packaging for more hygiene

EN 13034



Type PB [6]-B

EN 14126



EN ISO 13688

Industries

 Care Edition

 Medical Edition

Advantages

- + Breathable and water-repellent
- + Perfect hold
- + Protection against infection agents

Performance values

EN 13034:2005 + A1:2009	Test parameter	Performance levels	Test result
 Type PB [6]-B	Repellency against liquids: Sulphuric acid 30%	1-3	3
	Repellency against liquids: Sodium hydroxide 10%	1-3	3
	Repellency against liquids: o-Xylene	1-3	3
	Repellency against liquids: Butan-1-ol	1-3	3
	Resistance against penetration: Sulphuric acid 30%	1-3	3
	Resistance against penetration: Sodium hydroxide 10%	1-3	3
	Resistance against penetration: o-Xylene	1-3	3
	Resistance against penetration: Butan-1-ol	1-3	3
EN 14126:2003	Test parameter	Performance levels	Test result
	Resistance to penetration by contaminated liquids: synthetic blood	1-6	6
	Resistance to penetration by contaminated liquids: Phi-X174 bacteriophage	1-6	6
	Resistance to penetration by infective agents due to mechanical contact with substances containing contaminated liquids: Staphylococcus aureus	1-6	6
	Resistance to penetration by contaminated aerosols: Staphylococcus aureus	1-3	3
	Resistance to penetration by contaminated solid particles: Bacillus subtilis	1-3	3

Additional information

	Material		Colour		Grammage
Shell fabric	100% Polyethylene		light blue		
					
Do not wash	Do not bleach	Do not tumble dry	Do not iron	Do not dry clean	

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
2XL	4059904188273	4059904188280	0,5 m	0,32 m	0,35 m	13,88 kg	360

Customs tariff number: 62101098000

Country of origin: China