

007609

Consumables



Robust and durable
Easy to deform
Constant suction flow rate

007609

Material: Polypropylene, transparent

Miscellaneous: -

Packaging unit: 10 piece(s)

Saliva ejectors, transparent (colour code: 1400), polypropylene, robust and durable, length approx. 13 cm, removable cap, easy to deform, reliable suction of saliva and fluids, constant suction flow rate, suitable for standard extraction systems, polybag per 100 pieces

Industries

 Dental Edition

Advantages

- + Robust and durable
- + Easy to deform
- + Constant suction flow rate



Colour variants



007609
Colour code: 1400



007600
Colour code: 1100



007601
Colour code: 2200



007602
Colour code: 3000



007603
Colour code: 6200



007604
Colour code: 6300



007607
Colour code: 6000

Additional information

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904061378	4059904135413	0,3 m	0,19 m	0,22 m	4,24 kg	81

Customs tariff number: 39269097000

Country of origin: Germany