

143311 // FREEZE XL

Consumables
Medical device class IIa



Immediate and rapid cooling
Relieves pain and swelling
Versatile use at work and at home

143311 // FREEZE XL

Material: Plastic, white

Miscellaneous: -

Packaging unit: 20 piece(s)

NITRAS FREEZE XL, instant cold compress, white (colour code: 1100), Dimensions: 15 x 27 cm, for rapid cooling, ready for immediate use simply press to activate, relieves pain and swelling caused by sprains, strains, bruises, minor burns, headaches, toothaches and insect bites, flexible, supple shape, easy and quick to use, carton per 20 pieces

Industries

-  Beauty Edition
-  Black Edition
-  Care Edition
-  Dental Edition
-  Medical Edition

Advantages

- + Immediate and rapid cooling
- + Relieves pain and swelling
- + Versatile use at work and at home



Additional information**Logistics data**

Colour code	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
1100	4059904225756	4059904225145	0,34 m	0,27 m	0,17 m	5,5 kg	

Customs tariff number: 3824999699

Country of origin: China