

2008311 // BIG BLUE WAVE

Nitrile disposable gloves
PPE risk category III | Medical device class I





Very good tactile sensitivity
Very good dry and oil grip
Very robust and resistant



2008311 // BIG BLUE WAVE

Material:	Nitrile, blue
Coating:	-
Size(s):	S-XL
Packaging unit:	10 piece(s)
Certification:	EN ISO 21420, EN ISO 374, EN 421

NITRAS BIG BLUE WAVE, nitrile disposable gloves, blue (colour code: 2200), non-sterile, rolled edge, powder free, manufactured according to EN 455, for food contact, ambidextrous (left and right hand), medical gloves, micro-roughened fingertips, colour code system for sizes, AQL 1.5, protection against micro-organisms, bacteria and viruses, protection against radioactive contamination by solid particles, particularly sensitive, perfect for working with small parts, very good tactile sensitivity, very good dry and oil grip, very robust and resistant, medical devices class I, PPE risk category III, box per 200 pieces



Industries

-  Beauty Edition
-  Black Edition
-  Cleaning Edition
-  Dental Edition
-  Food Edition
-  Medical Edition

Advantages

- + Very good tactile sensitivity
- + Very good dry and oil grip
- + Very robust and resistant



Further images



Colour variants



2008311 // BIG BLUE WAVE
Colour code: 2200



2008310 // BIG WHITE WAVE
Colour code: 1100



2008320 // BIG BLACK WAVE
Colour code: 1000

Performance values

EN ISO 374-1:2016 + A1:2018	Test parameter	Performance levels	Test result
 JKT	J - n-Heptane	1-6	3 (> 60)
	K - Sodium hydroxide 40%	1-6	6 (> 480)
	T - Formaldehyde 37%	1-6	4 (> 120)
EN ISO 374-5:2016	Test parameter	Performance levels	Test result
 VIRUS	Penetration	-	-
	Protection against viruses	-	-
EN ISO 21420:2020 + A1:2024	Test parameter	Performance levels	Test result
EN 421:2010	Test parameter	Performance levels	Test result
			

Additional information

	Material	Colour
Palm	Nitrile	blue
Index finger and thumb	Nitrile	blue
Back of the hand	Nitrile	blue
Fingertips	Nitrile	blue
Cuff / knitted wrist	Nitrile	blue

Thickness (approx.)

Finger	0,10 mm
Palm	0,08 mm
Cuff	0,07 mm

Size	Length (mm)	Width palm (mm)
S	240	84
M	240	94
L	240	105
XL	240	113

Properties	Before ageing	After ageing
Elongation at break	500%	400%
Tear resistance	14 MPa	14 MPa

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
S	4059904037632	4059904006041	0,43 m	0,31 m	0,24 m	7,38 kg	60
M	4059904037649	4059904006058	0,43 m	0,31 m	0,24 m	8,09 kg	60
L	4059904037656	4059904006065	0,43 m	0,31 m	0,24 m	8,5 kg	60
XL	4059904037663	4059904006072	0,43 m	0,31 m	0,24 m	9,44 kg	60

Customs tariff number: 40151200000

Country of origin: Malaysia