

2008320 // BIG BLACK WAVE

Nitrile disposable gloves
PPE risk category III | Medical device class I



Very good tactile sensitivity
Very good dry and oil grip
Very robust and resistant



2008320 // BIG BLACK WAVE

Material:	Nitrile, black
Coating:	-
Size(s):	S-XL
Packaging unit:	10 piece(s)
Certification:	EN ISO 21420, EN ISO 374, EN 421

NITRAS BIG BLACK WAVE, nitrile disposable gloves, black (colour code: 1000), non-sterile, rolled edge, powder free, manufactured according to EN 455, for food contact, ambidextrous (left and right hand), medical gloves, micro-roughened fingertips, colour code system for sizes, AQL 1.5, protection against micro-organisms, bacteria and viruses, protection against radioactive contamination by solid particles, particularly sensitive, perfect for working with small parts, very good tactile sensitivity, very good dry and oil grip, very robust and resistant, medical devices class I, PPE risk category III, box per 200 pieces

ISO 374-1/Type B



JKT

ISO 374-5:2016



VIRUS

EN ISO 21420

EN 421:2010


AQL
1,5

NI

TESTED
ACCORDING TO
EN 455

Industries

-  Beauty Edition
-  Black Edition
-  Cleaning Edition
-  Dental Edition
-  Food Edition
-  Medical Edition

Advantages

- + Very good tactile sensitivity
- + Very good dry and oil grip
- + Very robust and resistant



Further images



Colour variants



2008320 // BIG BLACK WAVE
Colour code: 1000



2008310 // BIG WHITE WAVE
Colour code: 1100



2008311 // BIG BLUE WAVE
Colour code: 2200

Performance values

EN ISO 374-1:2016 + A1:2018	Test parameter	Performance levels	Test result
 JKT	J - n-Heptane	1-6	3 (> 60)
	K - Sodium hydroxide 40%	1-6	6 (> 480)
	T - Formaldehyde 37%	1-6	4 (> 120)
EN ISO 374-5:2016	Test parameter	Performance levels	Test result
 VIRUS	Penetration	-	Passed
	Protection against viruses	-	Passed
EN ISO 21420:2020 + A1:2024	Test parameter	Performance levels	Test result
	Dexterity	1-5	5
EN 421:2010	Test parameter	Performance levels	Test result
			

Additional information

	Material	Colour
Palm	Nitrile	black
Index finger and thumb	Nitrile	black
Back of the hand	Nitrile	black
Fingertips	Nitrile	black
Cuff / knitted wrist	Nitrile	black

Thickness (approx.)

Finger	0,10 mm
Palm	0,08 mm
Cuff	0,07 mm

Size	Length (mm)	Width palm (mm)
XS	240	76
S	240	84
M	240	94
L	240	105
XL	240	113

Properties	Before ageing	After ageing
Elongation at break	500%	400%
Tear resistance	14 MPa	14 MPa

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
XS	4059904126916	4059904126923	0,43	0,31	0,24 m	7,38 kg	60
S	4059904126930	4059904126947	0,43 m	0,31 m	0,24 m	7,38 kg	60
M	4059904126954	4059904126961	0,43 m	0,31 m	0,24 m	8,09 kg	60
L	4059904126978	4059904126985	0,43 m	0,31 m	0,24 m	8,5 kg	60
XL	4059904126992	4059904127005	0,43 m	0,31 m	0,24 m	9,44 kg	60

Customs tariff number: 40151200000

Country of origin: Malaysia