


20250

Sterilisation packaging



High tear and puncture resistance

Secure seal seams

Reliable sterile barrier

MADE IN
EUROPE

20250

Material: Polypropylene / Paper, blue

Miscellaneous: Width of roll 250 mm

Packaging unit: 2 piece(s)

Sterilisation foil with steam and EO gas indicator,, blue (colour code: 2200), polypropylene /paper, high tear and puncture resistance, roll length 200 m, roll width 250 mm, suitable for FO sterilisation, secure seal seams, reliable sterile barrier, sealing temperature: 160°C / 180°C, time between the machine's sealing bars and the product: 0.9 seconds / 1.5 seconds, pressure: 1–3 bar, recommended temperature: 160°C and 180°C, Steam sterilisation: 121°C = 15 minutes (minimum) / 134°C = 5 minutes (minimum), EtO gas: 50–60°C = 3–6 hours



Industries

-  Beauty Edition
-  Black Edition
-  Dental Edition
-  Medical Edition

Advantages

- + High tear and puncture resistance
- + Secure seal seams
- + Reliable sterile barrier



Additional information**Logistics data**

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904060982	4059904029798	0,43 m	0,23 m	0,29 m	11,54 kg	42

Customs tariff number: 48025525000

Country of origin: Italy