


210000

Sterilisation packaging



Plastic test specimens
Paper indicator strips
Use in autoclave

MADE IN
EUROPE

210000

Material: Plastics / Paper, white / blue

Miscellaneous: EN 867, EN ISO 11140

Packaging unit: 1 piece(s)

Helix-test for steam sterilisation, white / blue, (colour code: 1122), plastic test specimens, paper indicator strips, steam penetration test for hollow bodies for sterilisation process control, use in autoclave for 3.5 minutes at 134° C, set consisting of one test piece and 100 indicator strips per pack, complies with technical standards EN 867-5 and EN ISO 11140-1 Type 2

MADE IN
EUROPE

Industries

-  Beauty Edition
-  Black Edition
-  Dental Edition
-  Medical Edition

Advantages

- + Plastic test specimens
- + Paper indicator strips
- + Use in autoclave



Additional information

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size		4059904188266	0,16 m	0,09 m	0,03 m	0,2 kg	

Customs tariff number: 38221900190

Country of origin: Italy