

214260

Consumables



Sterilisation visible due to colour change
Can be used for several procedures
High level of safety for patients

MADE IN
EUROPE

214260

Material:	Polypropylene / Paper, blue
Miscellaneous:	Size 140 x 260 mm
Packaging unit:	8 piece(s)

Sterilisation bags with steam and EO gas indicator, blue (colour code: 2200), polypropylene / paper, high tear and puncture resistance, dimensions: 140 x 260 mm, self-adhesive, secure seal seams, reliable sterile barrier, box per 200 bags, the clear bags are used as packaging for medical products that are to be steam- or gas-sterilised. Seal the open end using a continuous or pulse sealer at the recommended temperature (170–180 °C) and pressure: 70–90 N/cm² (please check with the manufacturer of the sealing machine); the packaging must only be used once. Do not use damaged packaging. EO gas: yellow/brown (sterilised); steam: brown (sterilised). Avoid direct sunlight and excessive heat. Temperature: 10–30 °C; humidity: 30–60 °C

MADE IN
EUROPE

Industries

-  Beauty Edition
-  Black Edition
-  Dental Edition
-  Medical Edition

Advantages

- + Sterilisation visible due to colour change
- + Can be used for several procedures
- + High level of safety for patients



Additional information**Logistics data**

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904060999	4059904029804	0,34 m	0,34 m	0,26 m	8,6 kg	48

Customs tariff number: 48025525000

Country of origin: Italy