

# 275250

## Consumables



Sterilisation visible due to colour change  
Can be used for several procedures  
High level of safety for patients

MADE IN  
EUROPE

## 275250

|                        |                             |
|------------------------|-----------------------------|
| <b>Material:</b>       | Polypropylene / Paper, blue |
| <b>Miscellaneous:</b>  | Size 75 x 250 mm            |
| <b>Packaging unit:</b> | 12 piece(s)                 |

Sterilisation bags with steam and EO gas indicator, blue (colour code: 2200), polypropylene / paper, high tear and puncture resistance, dimensions: 75 x 250 mm, self-adhesive, secure seal seams, reliable sterile barrier, box per 200 bags, the clear bags are used as packaging for medical products that are to be steam- or gas-sterilised. Seal the open end using a continuous or pulse sealer at the recommended temperature (170–180 °C) and pressure: 70–90 N/cm<sup>2</sup> (please check with the manufacturer of the sealing machine); the packaging must only be used once. Do not use damaged packaging. EO gas: yellow/brown (sterilised); steam: brown (sterilised). Avoid direct sunlight and excessive heat. Temperature: 10–30 °C; humidity: 30–60 °C



### Industries

-  Beauty Edition
-  Black Edition
-  Dental Edition
-  Medical Edition

### Advantages

- + Sterilisation visible due to colour change
- + Can be used for several procedures
- + High level of safety for patients



## Additional information

### Logistics data

| Size     | GTIN (SPU)    | GTIN (PU)     | Length (m) | Width (m) | Height (m) | Weight (kg) | Quantity PU / pallet |
|----------|---------------|---------------|------------|-----------|------------|-------------|----------------------|
| One size | 4059904061019 | 4059904029828 | 0,34 m     | 0,33 m    | 0,26 m     | 7,18 kg     | 54                   |

Customs tariff number: 48025525000

Country of origin: Italy