

TECHNICAL DATASHEET

4015

Material: Polypropylene, white
Others: Diameter approx. 55 cm
Packaging unit: 1,000 pieces



Description

NITRAS clip caps, white, diameter approx. 55 cm, non-woven fabric, fibreglass-free and breathable, ruffled, good wearing comfort, double cotton-covered elastic band, 10 bags, 100 pieces each

Technical Specification

Weight: $12 \text{ g/sqm} \pm 2 \text{ g/sqm}$

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural materical changes during a period could implicate change in the properties.