

4227

Couch covers



Water-repellent
Tear-resistant and durable
Optimal hold

4227

Material: Chlorinated Polyethylene, pink

Miscellaneous: Size approx. 210 x 90 x 20 cm

Packaging unit: 10 piece(s)

NITRAS CPE mattress cover, pink (colour code: 6400), chlorinated polyethylene, tear-resistant and durable, water-repellent, thickness approx. 30 μ , dimensions: approx. 210 x 90 x 20 cm, with elastic band for optimal hold, easy to put on and take off, polybag per 10 pieces

Industries

 Beauty Edition

 Black Edition

Advantages

- + Water-repellent
- + Tear-resistant and durable
- + Optimal hold



Colour variants



4227
Colour code: 6400



4225
Colour code: 2200



4226
Colour code: 1000

Additional information**Logistics data**

Colour code	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
6400	4059904232822	4059904232839	0,5 m	0,3 m	0,42 m	6,04 kg	24

Customs tariff number: 39249000909

Country of origin: China