

TECHNICAL DATASHEET

4322 // SOFT PROTECT

Material:	Polypropylene, green	
Others:	One size	
Packaging Unit:	500 pieces	



Description

NITRAS SOFT PROTECT, medical face mask, made of very soft and fibreglass-free non-woven fabric, 3-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 Type IIR, bacterial filtration efficiency > 98%, 10 boxes, 50 pieces each

Construction

3-ply layer, made of polypropylene

Two elastic, latex-free elastic bands

Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

Technical specification

Layer 1: $23 \text{ g/sqm} \pm 2 \text{ g/sqm}$ Layer 2: $25 \text{ g/sqm} \pm 2 \text{ g/sqm}$ Layer 3: $23 \text{ g/sqm} \pm 2 \text{ g/sqm}$ Measurement: $17.5 \text{ cm} \times 9.5 \text{ cm} \pm 2 \text{ mm}$

Weight: $3 g \pm 0.2 g$

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural materical changes during a period could implicate change in the properties.