

4324 // SOFT PROTECT

Face masks
Medical device class I



Breathable whilst providing reliable protection
Integrated nosepiece for optimal fit
Latex-free ear loops for optimal hold

4324 // SOFT PROTECT

Material: Polypropylene, light pink

Miscellaneous: 3-ply

Packaging unit: 10 piece(s)

NITRAS SOFT PROTECT, medical face mask, light pink (colour code: 6300), non-woven fabric, 3-ply, fiberglass free, breathable whilst providing reliable protection, very soft, integrated nosepiece for optimal fit, round and latex-free ear loops for optimal hold, EN 14683 Type IIR, bacterial filtration efficiency > 98%, box per 50 pieces

EN 14683
Type IIR

Industries

-  Beauty Edition
-  Care Edition
-  Dental Edition
-  Medical Edition

Advantages

- + Breathable whilst providing reliable protection
- + Integrated nosepiece for optimal fit
- + Latex-free ear loops for optimal hold



Further images



Colour variants



4324 // SOFT PROTECT
Colour code: 6300



4320 // SOFT PROTECT
Colour code: 1100



4321 // SOFT PROTECT
Colour code: 2300



4322 // SOFT PROTECT
Colour code: 3000



4323 // SOFT PROTECT
Colour code: 6500



4329 // SOFT PROTECT
Colour code: 6400

Additional information

	Material	Colour	Grammage
Shell fabric	100% Polypropylene	light pink	

Physical properties

Grammage	23 g/qm
Dimensions	17,5 cm x 9,5 cm ± 2 mm

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904073340	4059904073357	0,53 m	0,2 m	0,21 m	2,42 kg	72

Customs tariff number: 63079095110

Country of origin: China