

TECHNICAL DATASHEET

432800 // SOFT PROTECT VOGUE

Material: Polypropylen, white / black
Others: One size
Packaging Unit: 500 pieces



Description

NITRAS MEDICAL PROTECT, medical face mask, made of fibreglass-free non-woven fabric, 3-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 type IIR, bacterial filtration efficiency > 98%, 10 boxes, 50 pieces each

Construction

Layer 1: patterned colored polypropylene layer

Layer 2: white Meltblown-filter

Layer 3: white layer made of polypropylene

Two elastic, latex-free elastic bands

Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

Technical specification

 $\begin{array}{lll} \mbox{Layer 1:} & 23 \mbox{ g/sqm} \pm 2 \mbox{ g/sqm} \\ \mbox{Layer 2:} & 25 \mbox{ g/sqm} \pm 2 \mbox{ g/sqm} \\ \mbox{Layer 3:} & 23 \mbox{ g/sqm} \pm 2 \mbox{ g/sqm} \\ \mbox{Measurement:} & 17.5 \mbox{ cm} \times 9.5 \mbox{ cm} \pm 2 \mbox{ mm} \\ \end{array}$

Weight: $3.0 \text{ g} \pm 0.2 \text{ g}$

The medical face mask offers good protection against dust, splash water and bacteria and is specially designed for medical use.

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural materical changes during a period could implicate change in the properties.