

TECHNICAL DATASHEET

4331HY // PROTECT

Material:	Polypropylene, blue
Others:	One size
Packaging Unit:	2,000 pieces



Description

NITRAS PROTECT, medical face mask, 3-ply, integrated nosepiece, round and latex-free elastic bands, EN 14683 type IIR, , bacterial filtration efficiency > 98%, 40 boxes, 50 pieces each

Construction

3-ply blue layer, made of polypropylene
Two elastic, latex-free elastic bands
Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

Technical specification

Layer 1: 15 g/sqm \pm 2 g/sqm
Layer 2: 25 g/sqm \pm 2 g/sqm
Layer 3: 30 g/sqm \pm 2 g/sqm
Measurement: 17.5 cm x 9.5 cm \pm 2 mm
Weight: 2.3 g \pm 0.2 g

The medical face mask offers good protection against dust, splash water and bacteria and was specially designed for medical use.

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural material changes during a period could implicate change in the properties.