

## TECHNICAL DATASHEET

### 4331 // PROTECT

Material:	Polypropylene, blue
Others:	One size
Packaging Unit:	500 pieces



### Description

NITRAS MEDICAL PROTECT, medical face mask, made of fibreglass-free non-woven fabric, 3-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 type IIR, bacterial filtration efficiency > 98%, 10 boxes, 50 pieces each

### Construction

3-ply blue layer, made of polypropylene  
 Two elastic, latex-free elastic bands  
 Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

### Technical specification

Layer 1:	15 g/sqm ± 2 g/sqm
Layer 2:	25 g/sqm ± 2 g/sqm
Layer 3:	30 g/sqm ± 2 g/sqm
Measurement:	17.5 cm x 9.5 cm ± 2 mm
Weight :	2.3 g ± 0.2 g

### Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural material changes during a period could implicate change in the properties.