

TECHNICAL DATASHEET

4332 // PROTECT

Material: Others: Packaging Unit: Polypropylene, green One size 500 pieces



Description

NITRAS MEDICAL PROTECT, medical face mask, green, made of fibreglass-free non-woven fabric, 3-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 type IIR, 10 boxes, 50 pieces each

Construction

3-ply green layer, made of polypropylene Two elastic, latex-free elastic bands Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

Technical specification

Layer 1:	15 g/sqm \pm 2 g/sqm
Layer 2:	25 g/sqm \pm 2 g/sqm
Layer 3:	30 g/sqm \pm 2 g/sqm
Measurement:	17.5 cm x 9.5 cm \pm 2 mm
Weight :	$2.3~{ m g}\pm 0.2~{ m g}$

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural materical changes during a period could implicate change in the properties.