

TECHNICAL DATASHEET

4338 // PROTECT

| Material: | Polypropylene, black | |
|-----------------|----------------------|--|
| Others: | One size | |
| Packaging Unit: | 500 pieces | |



Description

NITRAS MEDICAL PROTECT, medical face mask, black, made of fibreglass-free non-woven fabric, 3-ply, integrated nosepiece, round and latex-free elastic bands, manufactured according to EN 14683 type IIR, 10 boxes, 50 pieces each

Construction

3-ply black layer, made of polypropylene

Two elastic, latex-free elastic bands

Nosepiece, made of plastic coated wire, that can be adjusted on the bridge of the nose

Technical specification

 $\begin{array}{lll} \mbox{Layer 1:} & \mbox{15 g/sqm} \pm 2 \mbox{g/sqm} \\ \mbox{Layer 2:} & \mbox{25 g/sqm} \pm 2 \mbox{g/sqm} \\ \mbox{Layer 3:} & \mbox{30 g/sqm} \pm 2 \mbox{g/sqm} \\ \mbox{Measurement:} & \mbox{17.5 cm} \times 9.5 \mbox{cm} \pm 2 \mbox{mm} \\ \end{array}$

Weight: $2.3 g \pm 0.2 g$

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural materical changes during a period could implicate change in the properties.