

500057

Accessories



Very robust
Reliable protection against drying out
Refillable and reusable

500057

Material: Polyethylene, white / blue

Miscellaneous: -

Packaging unit: 1 piece(s)

Dispenser bucket for wet non-woven disinfection wipes, blue / white (colour code: 2211), polyethylene, very robust, with carrying handle, transparent lid with blue dispensing opening, reliable protection against drying out, capacity approx. 5,700 ml, height approx. 19.4 cm, refillable and reusable, easy and hygienic dispensing of wipes

Industries

-  Care Edition
-  Dental Edition
-  Medical Edition

Advantages

- + Very robust
- + Reliable protection against drying out
- + Refillable and reusable



Additional information

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size		4059904125599	0,24 m	0,23 m	0,22 m	0,25 kg	1080

Customs tariff number: 39269097900

Country of origin: Germany