

51552

Consumables



Puncture-resistant and sturdy

Secure-locking lid

Capacity 2,000 ml

51552

Material: Polyethylene, black

Miscellaneous: EN ISO 23907

Packaging unit: 30 piece(s)

NITRAS sharps container, black (colour code: 1000), polyethylene, puncture-resistant and sturdy, dimensions: approx. 20 x 16 x 12 cm, capacity 2,000 ml, secure-locking lid, transparent with red opening, EN ISO 23907

Industries

-  Black Edition
-  Medical Edition

Advantages

- + Puncture-resistant and sturdy
- + Secure-locking lid
- + Capacity 2,000 ml



Additional information**Logistics data**

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904173132	4059904173149	0,41 m	0,34 m	0,25 m	5,5 kg	48

Customs tariff number: 39231090900

Country of origin: China