

51655

Consumables



Puncture-resistant and sturdy
Secure-locking lid
Capacity 5,000 ml

51655

Material: Polyethylene, black

Miscellaneous: EN ISO 23907

Packaging unit: 24 piece(s)

NITRAS sharps container, black (colour code: 1000), polyethylene, puncture-resistant and sturdy, dimensions: approx. 24.5 x 20 x 16.5 cm, capacity 5,000 ml, secure-locking lid, transparent with red opening, EN ISO 23907

Industries

-  Black Edition
-  Medical Edition

Advantages

- + Puncture-resistant and sturdy
- + Secure-locking lid
- + Capacity 5,000 ml



Additional information

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904098466	4059904092433	0,52 m	0,44 m	0,3 m	6,64 kg	24

Customs tariff number: 39231090900

Country of origin: China