

# 51655

Consumables



Puncture-resistant and sturdy  
Secure-locking lid  
Capacity 5,000 ml

## 51655

**Material:** Polyethylene, black

**Miscellaneous:** EN ISO 23907

**Packaging unit:** 24 piece(s)

NITRAS sharps container, black (colour code: 1000), polyethylene, puncture-resistant and sturdy, dimensions: approx. 24.5 x 20 x 16.5 cm, capacity 5,000 ml, secure-locking lid, transparent with red opening, EN ISO 23907

### Industries

-  Black Edition
-  Medical Edition

### Advantages

- + Puncture-resistant and sturdy
- + Secure-locking lid
- + Capacity 5,000 ml



**Additional information****Logistics data**

| Size     | GTIN (SPU)    | GTIN (PU)     | Length (m) | Width (m) | Height (m) | Weight (kg) | Quantity PU /<br>pallet |
|----------|---------------|---------------|------------|-----------|------------|-------------|-------------------------|
| One size | 4059904098466 | 4059904092433 | 0,52 m     | 0,44 m    | 0,3 m      | 6,64 kg     | 24                      |

Customs tariff number: 39231090900

Country of origin: China