

56070

Accessories



Lint-free
Antiseptic
Fast and reliable disinfection

56070

Material: Viscose / polyethylenterephthalate, white

Miscellaneous: -

Packaging unit: 20 piece(s)

Alcohol swabs, white (colour code: 1100), viscose / polyethylenterephthalate, soaked with 70% isopropyl alcohol, lint-free, antiseptic, dimensions: approx. 6.5 x 3 cm, folded approx. 3 x 3 cm, ready-to-use, fast and reliable disinfection, single packaging for more hygiene, box per 100 pieces

Industries

-  Beauty Edition
-  Black Edition
-  Medical Edition

Advantages

- + Lint-free
- + Antiseptic
- + Fast and reliable disinfection



Additional information**Logistics data**

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904064904	4059904064911	0,28 m	0,22 m	0,13 m	2,14 kg	240

Customs tariff number: 63071030000

Country of origin: China