

# TECHNICAL DATA SHEET

# 58100 | 58112 // PROTECTASEPT

Туре:	Disinfection wipes, ready-to-use, Lemon	
Dimensions:	13 x 20 cm	
Packaging unit:	1.200 pieces (12 cans / refill bags 100 pieces each)	



ARTICLE	BUNDLES	PACKAGING UNIT
58100	Can	12 cans
58112	Refill bags	12 refill bags

#### Description

Ready-to-use disinfectant wipes for universal disinfection of alcohol-compatible surfaces, fragrance: lemon, IHO-listed, wipe size 12 x 20 cm, box / bag of 100 wipes made of PP/VIS mixed material, impregnated with wipe disinfectant, ready-to-use, efficacy according to EN 16615

## Properties

Ready-to-use alcohol wipes - selected lint-free nonwoven qualities for very good microbiological values in the EN16615 test method with consolidated fibres, extremely broad spectrum of efficacy for the highest requirements within 2 minutes at high loads, practical and easy application and removal from the dispenser box, "One-hand removal" due to special sheet perforation with large and pre-plucked roll core, practical printing for maximum regulatory safety with minimum documentation, long service life in dispenser box (assessed), evaporates quickly and leaves no residues, approval as biocide and medical device (dual use)

#### **Preparation type**

Ready-to-use, pre-soaked wipes for the rapid disinfection of medical inventory as well as near-patient and hygienically critical areas. In the practical dispenser box, they can be removed individually and are quickly ready for use.

#### **Areas of application**

According to MPG: Alcoholic wet wipes for prophylactic surface disinfection of alcohol-resistant medical devices. Not suitable for final (terminal) disinfection of invasive medical devices.

According to the Biocidal Products Ordinance: Alcoholic wet wipes for rapid disinfection of alcohol-resistant surfaces.

#### Composition

In 100g impregnating liquid are contained: 43.0g ethanol; 8.0g 1-propanol; 0.055g didecyldimethyl-ammonium chloride.

 Kreissparkasse Köln

 IBAN:
 DE49 3705 0299 0127 2721 64

 BIC:
 COKSDE33XXX

 Sparkasse Köln/Bonn

 IBAN:
 DE32 3705 0198 1933 3245 90

 BIC:
 COLSDE33XXX

 sc de odge auf Anfrage)
 Our general term

 HypoVereinsbank

 IBAN:
 DE63 3702 0090 0027 7497 39

 BIC:
 HYVEDEMM429

Es gelten unsere allgemeinen Geschäfts- und Lieferbedingungen (zu finden unter www.nitras.de oder auf Anfrage). Our general terms and conditions of business and delivery apply (can be found at www.nitras.de or on request).

MCD Medical Care Dental GmbH | Heinrich-Hertz-Str. 7 | 50181 Bedburg | Germany | Geschäftsführer / Managing Director: Maria Baric, Markus Seiler | HRB 83971 | Gerichtsstand Köln | USt.-IdNr.: DE299692222

🌭 +49 2272 401 111

www.nitras-medical.com



### Application

Wipe surfaces and objects until they are completely wetted and allow the disinfectant to act. After the exposure time has elapsed, wipe with a disposable cloth if necessary. Only suitable for use on alcohol-compatible materials (unsuitable for acrylic glass). The microbiological requirements for the quality of the rinsing water are at least drinking water quality.

#### **Exposure time**

Contact times for cloth wiping systems: Bacteria and Candida albicans according to EN 16615 with high load: Conc. 2 min.

Relative to the impregnating liquid: Bacteria according to EN 13727 and EN 13697 (high load): Conc. 1 min. Candida albicans according to EN 13624 and EN 13697 (high load): Conc. 1min. TbB inactivation according to EN 14348 (M.terrae) with high load: conc. 1min. Enveloped viruses according to DVV/2012 in the area test with high load: Conc. 1min. Adenoviruses according to EN 14476 with high load: Conc. 1min. Noro viruses according to EN 14476 with high load: conc. 1min. Rota viruses according to EN 14476 with high load: Conc. 1min.

#### Additional information on the spectrum of activity

Bactericidal (incl. MRSA), levurocidal (Candida albicans), tuberculocidal (M. terrae), "limited virucidal PLUS" (effective against enveloped viruses such as e.g. HBV, HIV, HCV) and non-enveloped adenoviruses, noroviruses, rotaviruses

#### **Safety instructions**

H 228 Flammable solid. Keep away from heat / sparks / open flame / hot surfaces. Do not smoke. Wear protective gloves / eye protection. Keep in a closed container.



# Distribution

Distribution according to the Biocidal Products Directive only in Germany and Slovenia

Kreissparkasse Köln IBAN: DE49 3705 0299 0127 2721 64 BIC: COKSDE33XXX

Sparkasse Köln/Bonn IBAN: DE32 3705 0198 1933 3245 90 BIC: COLSDE33XXX

HypoVereinsbank IBAN: DE63 3702 0090 0027 7497 39 BIC: HYVEDEMM429

Es gelten unsere allgemeinen Geschäfts- und Lieferbedingungen (zu finden unter www.nitras.de oder auf Anfrage). Our general terms and conditions of business and delivery apply (can be found at www.nitras.de or on request).

MCD Medical Care Dental GmbH | Heinrich-Hertz-Str. 7 | 50181 Bedburg | Germany | Geschäftsführer / Managing Director: Maria Baric, Markus Seiler | HRB 83971 | Gerichtsstand Köln | USt.-IdNr.: DE299692222

+49 2272 401 111

www.nitras-medical.com