

76050

Couch covers



Absorbent top layer
Liquid-impermeable underside
Protection against humidity and dirt

76050

Material: Paper / Polyethylene, black

Miscellaneous: 3-ply

Packaging unit: 6 piece(s)

Couch cover, black (colour code: 1000), paper / polyethylene, tear-resistant and durable, 3-ply: two layers paper, one layer polyethylene, absorbent top layer, liquid-impermeable underside, perforated, sheet width approx. 60 cm, sheet length approx. 50 cm, protection against humidity and dirt, roll length 40 m

Industries

 Black Edition

Advantages

- + Absorbent top layer
- + Liquid-impermeable underside
- + Protection against humidity and dirt



Additional information**Logistics data**

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size	4059904093737	4059904093744	0,63 m	0,34 m	0,22 m	7,82 kg	28

Customs tariff number: 48189010

Country of origin: Germany