

772501

Cups



Supports independent fluid intake
Reduces spillage of drinks
Lightweight and shatterproof

772501

Material: Polypropylene, blue

Miscellaneous: -

Packaging unit: 1 piece(s)

NITRAS beak cup, blue (colour code: 2200), polypropylene, robust and durable, 3-piece, capacity 250 ml, one lid with 12x10 mm opening, one lid with 4x4 mm opening, drinking spout for controlled drinking, supports independent fluid intake, reduces spillage of drinks, lightweight and shatterproof, easy to grip and handle, easy to clean, tested according to EN 12875-1 and 12875-2, carton per 50 sets



MADE IN
EUROPE

Industries



Care Edition

Advantages

- + Supports independent fluid intake
- + Reduces spillage of drinks
- + Lightweight and shatterproof



Colour variants



772501
Colour code: 2200



772500
Colour code: 1400



772502
Colour code: 3000



772505
Colour code: 4300



772507
Colour code: 6000

Additional information**Logistics data**

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
One size		4059904188044	0,39 m	0,27 m	0,14 m	1,75 kg	112

Customs tariff number: 39231090000

Country of origin: Croatia