

772502

Cups



Supports independent fluid intake
Reduces spillage of drinks
Lightweight and shatterproof



772502

Material: Polypropylene, green

Miscellaneous: -

Packaging unit: 1 piece(s)

NITRAS beak cup, green (colour code: 3000), polypropylene, robust and durable, 3-piece, capacity 250 ml, one lid with 12x10 mm opening, one lid with 4x4 mm opening, drinking spout for controlled drinking, supports independent fluid intake, reduces spillage of drinks, lightweight and shatterproof, easy to grip and handle, easy to clean, tested according to EN 12875-1 and 12875-2, carton per 50 sets



MADE IN EUROPE

Industries

 Care Edition

Advantages

- + Supports independent fluid intake
- + Reduces spillage of drinks
- + Lightweight and shatterproof



Colour variants



772502
Colour code: 3000



772500
Colour code: 1400



772501
Colour code: 2200



772505
Colour code: 4300



772507
Colour code: 6000

Additional information

Logistics data

| Size | GTIN (SPU) | GTIN (PU) | Length (m) | Width (m) | Height (m) | Weight (kg) | Quantity PU / pallet |
|----------|------------|---------------|------------|-----------|------------|-------------|----------------------|
| One size | | 4059904188051 | 0,39 m | 0,27 m | 0,14 m | 1,75 kg | 112 |

Customs tariff number: 39231090000

Country of origin: Croatia