

8100 // VINYL PF

Vinyl disposable gloves
PPE risk category I | Medical device class I




Particularly sensitive
Perfect for working with small parts
Very good tactile sensitivity



8100 // VINYL PF

Material:	Vinyl, transparent
Coating:	-
Size(s):	XS-XL
Packaging unit:	10 piece(s)
Certification:	-

NITRAS VINYL PF, vinyl disposable gloves, transparent (colour code: 1400), non-sterile, rolled edge, powder free, manufactured according to EN 455, for food contact, ambidextrous (left and right hand), medical gloves, AQL 1.5, particularly sensitive, perfect for working with small parts, very good tactile sensitivity, medical devices class I, PPE risk category I, box per 100 pieces



Industries

-  Care Edition
-  Cleaning Edition
-  Food Edition
-  Medical Edition

Advantages

- + Particularly sensitive
- + Perfect for working with small parts
- + Very good tactile sensitivity



Further images



Additional information

	Material	Colour
Palm	Vinyl	transparent
Index finger and thumb	Vinyl	transparent
Back of the hand	Vinyl	transparent
Fingertips	Vinyl	transparent
Cuff / knitted wrist	Vinyl	transparent

Thickness (approx.)

Finger	0,06 mm
Palm	0,06 mm
Cuff	0,06 mm

Size	Length (mm)	Width palm (mm)
XS	240	75
S	240	85
M	240	95
L	240	105
XL	240	115

Properties	Before ageing	After ageing
Elongation at break	300%	
Tear resistance	11 MPa	

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
XS	4059904092129	4059904092136	0,35 m	0,27 m	0,24 m	4,9 kg	80
S	4059904092082	4059904092099	0,35	0,27	0,24 m	5,49 kg	80
M	4059904092068	4059904092075	0,35	0,27	0,24 m	5,76 kg	80
L	4059904092044	4059904092051	0,35	0,27	0,24 m	6,72 kg	80
XL	4059904092105	4059904092112	0,35 m	0,27 m	0,24 m	6,48 kg	80

Customs tariff number: 39262000000

Country of origin: China