

8101 // VINYL P

Vinyl disposable gloves
PPE risk category I | Medical device class I




Particularly sensitive
Perfect for working with small parts
Very good tactile sensitivity



8101 // VINYL P

Material:	Vinyl, transparent
Coating:	-
Size(s):	S-XL
Packaging unit:	10 piece(s)
Certification:	-

NITRAS VINYL P, vinyl disposable gloves, transparent (colour code: 1400), non-sterile, rolled edge, powdered, manufactured according to EN 455, for food contact, ambidextrous (left and right hand), medical gloves, AQL 1.5, particularly sensitive, perfect for working with small parts, very good tactile sensitivity, medical devices class I, PPE risk category I, box per 100 pieces



Industries

-  Care Edition
-  Cleaning Edition
-  Food Edition
-  Medical Edition

Advantages

- + Particularly sensitive
- + Perfect for working with small parts
- + Very good tactile sensitivity



Further images



Additional information

	Material	Colour
Palm	Vinyl	transparent
Index finger and thumb	Vinyl	transparent
Back of the hand	Vinyl	transparent
Fingertips	Vinyl	transparent
Cuff / knitted wrist	Vinyl	transparent

Thickness (approx.)

Finger	0,6 mm
Palm	0,6 mm
Cuff	0,6 mm

Size	Length (mm)	Width palm (mm)
S	240	85
M	240	95
L	240	105
XL	240	115

Properties	Before ageing	After ageing
Elongation at break	300%	
Tear resistance	11 MPa	

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
S	4059904092181	4059904092198	0,35 m	0,27 m	0,24 m	5,08 kg	80
M	4059904092167	4059904092174	0,35 m	0,27 m	0,24 m	6,42 kg	80
L	4059904092143	4059904092150	0,35 m	0,27 m	0,24 m	6,26 kg	80
XL	4059904092204	4059904092211	0,35 m	0,27 m	0,24 m	6,73 kg	80

Customs tariff number: 39262000000

Country of origin: China