

TECHNICAL DATASHEET

8101 // VINYL P

Material:	Vinyl, transparent
Coating:	-
Size(s):	XS-XL
Packaging unit:	1,000 pieces
Certification:	-



Description

NITRAS VINYL P, vinyl disposable gloves, transparent, non-sterile, rolled edge, powdered, manufactured according to EN 455, for food contact, ambidextrous, medical examination gloves, AQL 1.5, 10 boxes, 100 pieces each

	Material	Colour
Palm	Vinyl	Transparent
Index finger and thumb	Vinyl	Transparent
Back of the hand	Vinyl	Transparent
Fingertips	Vinyl	Transparent
Cuff / knitted wrist	Vinyl	Transparent

Thickness (approx.)	mm
Finger	0,06
Palm	0,06
Cuff	0,06

Dimensions	Size	Values
Length (mm)	all sizes	min. 240
	XS	75 ± 5
Width of palm (mm)	S	85 ± 5
	M	95 ± 5
	L	105 ± 5
	XL	115 ± 5

Material properties	Standard before ageing
Elongation at break (%)	min. 300
Tensile strength (MPa)	min. 11

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural material changes during a period could implicate change in the properties.

Durability

See overpack