

TECHNICAL DATASHEET

8210 // SUPREME LATEX

Material:	Latex , natural colour
Coating:	-
Size(s):	XS-XL
Packaging unit:	1,000 pieces
Certification:	-



Description

NITRAS SUPREME LATEX, latex disposable gloves, natural colour, non-sterile, rolled edge, powder-free, manufactured according to EN 455, for food contact, ambidextrous, medical examination gloves, double chlorinated, micro-roughened surface, colour code system for sizes, AQL 1.5, 10 boxes, 100 pieces each

	Material	Colour
Palm	Latex	Natural colour
Index finger and thumb	Latex	Natural colour
Back of the hand	Latex	Natural colour
Fingertips	Latex	Natural colour
Cuff / knitted wrist	Latex	Natural colour

Thickness (approx.)	mm
Finger	0,17 ± 0,03
Palm	0,14 ± 0,02
Cuff	0,10 ± 0,02

Dimensions	Size	Values
Length (mm)	all sizes	min. 230
	XS	78 ± 3
Width of palm (mm)	S	83 ± 3
	M	93 ± 3
	L	102 ± 3
	XL	112 ± 3

Material properties	Standard before ageing	Standard after ageing
Elongation at break (%)	min. 750 (Guideline: 650-750)	min. 650 (Guideline: 500-600)
Tensile strength (MPa)	min. 21 (Guideline: 18-22)	min. 18 (Guideline: 14-18)

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural material changes during a period could implicate change in the properties.

Durability

See overpack