

# 8215 // POLYMER SOFT

Latex disposable gloves  
PPE risk category I | Medical device class I



Micro-roughened surface  
Particularly sensitive  
Very good tactile sensitivity



AQL  
1,5

LA

TESTED  
ACCORDING TO  
EN 455

## 8215 // POLYMER SOFT

|                 |              |
|-----------------|--------------|
| Material:       | Latex, white |
| Coating:        | -            |
| Size(s):        | XS-XL        |
| Packaging unit: | 10 piece(s)  |
| Certification:  | -            |

NITRAS POLYMER SOFT, latex disposable gloves, white (colour code: 1100), non-sterile, rolled edge, powder free, manufactured according to EN 455, extra grip, ambidextrous (left and right hand), micro-roughened surface, polymer coating, medical gloves, AQL 1.5, particularly sensitive, perfect for working with small parts, very slip resistant, very good tactile sensitivity, very good dry and wet grip, medical devices class I, PPE risk category I, box per 100 pieces



### Industries

-  Dental Edition
-  Medical Edition

### Advantages

- + Micro-roughened surface
- + Particularly sensitive
- + Very good tactile sensitivity



## Further images

---



## Additional information

|                        | Material | Colour |
|------------------------|----------|--------|
| Palm                   | Latex    | white  |
| Index finger and thumb | Latex    | white  |
| Back of the hand       | Latex    | white  |
| Fingertips             | Latex    | white  |
| Cuff / knitted wrist   | Latex    | white  |

### Thickness (approx.)

|        |       |
|--------|-------|
| Finger | 12 mm |
| Palm   | 11 mm |
| Cuff   | 11 mm |

| Size | Length (mm) | Width palm (mm) |
|------|-------------|-----------------|
| XS   | 240         | 76              |
| S    | 240         | 84              |
| M    | 240         | 94              |
| L    | 240         | 105             |
| XL   | 240         | 113             |

| Properties          | Before ageing | After ageing |
|---------------------|---------------|--------------|
| Elongation at break | 650%          | 500%         |
| Tear resistance     | 18 MPa        | 14 MPa       |

## Logistics data

| Size | GTIN (SPU)    | GTIN (PU)     | Length (m) | Width (m) | Height (m) | Weight (kg) | Quantity PU / pallet |
|------|---------------|---------------|------------|-----------|------------|-------------|----------------------|
| XS   | 4059904089730 | 4059904089785 | 0,37 m     | 0,25 m    | 0,23 m     | 5,83 kg     | 80                   |
| S    | 4059904089747 | 4059904089792 | 0,37 m     | 0,25 m    | 0,23 m     | 6,11 kg     | 80                   |
| M    | 4059904089754 | 4059904089808 | 0,37 m     | 0,25 m    | 0,23 m     | 6,76 kg     | 80                   |
| L    | 4059904089761 | 4059904089815 | 0,37 m     | 0,25 m    | 0,23 m     | 7,1 kg      | 80                   |
| XL   | 4059904089778 | 4059904089822 | 0,37 m     | 0,25 m    | 0,23 m     | 7,17 kg     | 80                   |

Customs tariff number: 40151200000

Country of origin: Malaysia