

8215 // POLYMER SOFT

Latex disposable gloves
PPE risk category I | Medical device class I



Micro-roughened surface
Particularly sensitive
Very good tactile sensitivity



AQL
1,5

LA

TESTED
ACCORDING TO
EN 455

8215 // POLYMER SOFT

Material:	Latex, white
Coating:	-
Size(s):	XS-XL
Packaging unit:	10 piece(s)
Certification:	-

NITRAS POLYMER SOFT, latex disposable gloves, white (colour code: 1100), non-sterile, rolled edge, powder free, manufactured according to EN 455, extra grip, ambidextrous (left and right hand), micro-roughened surface, polymer coating, medical gloves, AQL 1.5, particularly sensitive, perfect for working with small parts, very slip resistant, very good tactile sensitivity, very good dry and wet grip, medical devices class I, PPE risk category I, box per 100 pieces



Industries

-  Dental Edition
-  Medical Edition

Advantages

- + Micro-roughened surface
- + Particularly sensitive
- + Very good tactile sensitivity



Further images



Additional information

	Material	Colour
Palm	Latex	white
Index finger and thumb	Latex	white
Back of the hand	Latex	white
Fingertips	Latex	white
Cuff / knitted wrist	Latex	white

Thickness (approx.)

Finger	12 mm
Palm	11 mm
Cuff	11 mm

Size	Length (mm)	Width palm (mm)
XS	240	76
S	240	84
M	240	94
L	240	105
XL	240	113

Properties	Before ageing	After ageing
Elongation at break	650%	500%
Tear resistance	18 MPa	14 MPa

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
XS	4059904089730	4059904089785	0,37 m	0,25 m	0,23 m	5,83 kg	80
S	4059904089747	4059904089792	0,37 m	0,25 m	0,23 m	6,11 kg	80
M	4059904089754	4059904089808	0,37 m	0,25 m	0,23 m	6,76 kg	80
L	4059904089761	4059904089815	0,37 m	0,25 m	0,23 m	7,1 kg	80
XL	4059904089778	4059904089822	0,37 m	0,25 m	0,23 m	7,17 kg	80

Customs tariff number: 40151200000

Country of origin: Malaysia