

TECHNICAL DATASHEET

8215 // POLYMER SOFT

Material:	Latex , white
Coating:	-
Size(s):	XS-XL
Packaging unit:	1,000 pieces
Certification:	-



Description

NITRAS POLYMER SOFT, latex disposable gloves, white, non-sterile, rolled edge, powder-free, manufactured according to EN 455, ambidextrous, medical examination gloves, perfect grip, polymer coating, micro-roughened surface, AQL 1.5, 10 boxes, 100 pieces each

	Material	Colour
Palm	Latex	White
Index finger and thumb	Latex	White
Back of the hand	Latex	White
Fingertips	Latex	White
Cuff / knitted wrist	Latex	White

Thickness (approx.)	mm
Finger	0,12 ± 0,02
Palm	0,11 ± 0,02

Dimensions	Size	Values
Length (mm)	all sizes	min. 240
	XS	76 ± 3
Width of palm (mm)	S	84 ± 3
	M	94 ± 3
	L	105 ± 3
	XL	113 ± 3

Material properties	Standard before ageing	Standard after ageing
Elongation at break (%)	min. 650	min. 500
Tensile strength (MPa)	min. 18	min. 14

Storage

The products must be stored appropriately, this means in cartons in cool, dry rooms. External influences such as moisture, temperatures, light, as well as natural material changes during a period could implicate change in the properties.

Durability

See overpack