

8220 // BLACK SCORPION

Latex disposable gloves
PPE risk category III | Medical device class I



Very slip resistant
Very good tactile sensitivity
Very good dry and wet grip



8220 // BLACK SCORPION

Material:	Latex, black
Coating:	-
Size(s):	XS-XL
Packaging unit:	10 piece(s)
Certification:	EN ISO 21420, EN ISO 374

NITRAS BLACK SCORPION, latex disposable gloves, black (colour code: 1000), non-sterile, rolled edge, powder free, manufactured according to EN 455, for food contact, ambidextrous (left and right hand), medical gloves, extra grip, micro-roughened surface, colour code system for sizes, AQL 1.5, protection against micro-organisms and bacteria, particularly sensitive, perfect for working with small parts, very slip resistant, very good tactile sensitivity, very good dry and wet grip, medical devices class I, PPE risk category III, box per 100 pieces

ISO 374-1/Type C



ISO 374-5:2016



EN ISO 21420

LA


AQL
1,5

TESTED
ACCORDING TO
EN 455

Industries

-  Black Edition
-  Food Edition

Advantages

- + Very slip resistant
- + Very good tactile sensitivity
- + Very good dry and wet grip



Further images



Performance values

EN ISO 374-1:2016 + A1:2018	Test parameter	Performance levels	Test result
	K - Sodium hydroxide 40%	1-6	6 (> 480)
	L - Sulphuric acid 96%	1-6	1 (> 10)
	T - Formaldehyde 37%	1-6	6 (> 480)
EN ISO 374-5:2016	Test parameter	Performance levels	Test result
	Penetration	-	Passed
	Protection against viruses	-	-
EN ISO 21420:2020	Test parameter	Performance levels	Test result
	Dexterity	1-5	5

Additional information

	Material	Colour
Palm	Latex	black
Index finger and thumb	Latex	black
Back of the hand	Latex	black
Fingertips	Latex	black
Cuff / knitted wrist	Latex	black

Thickness (approx.)

Finger	0,17 mm
Palm	0,14 mm
Cuff	0,10 mm

Size	Length (mm)	Width palm (mm)
XS	240	76
S	240	84
M	240	94
L	240	105
XL	240	113

Properties	Before ageing	After ageing
Elongation at break	750%	240%
Tear resistance	21 MPa	18 MPa

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
XS	4059904086500	4059904086555	0,41 m	0,27 m	0,25 m	6,1 kg	72
S	4059904086517	4059904086562	0,4 m	0,27 m	0,25 m	6,81 kg	72
M	4059904086524	4059904086579	0,4 m	0,27 m	0,25 m	7,23 kg	72
L	4059904086531	4059904086586	0,4 m	0,27 m	0,25 m	8,56 kg	72
XL	4059904086548	4059904086593	0,4 m	0,27 m	0,25 m	8,96 kg	72

Customs tariff number: 40151200000

Country of origin: Malaysia