

8305 // PROTECT 300

Nitrile disposable gloves
PPE risk category III | Medical device class I



Very good tactile sensitivity
Very good dry and oil grip
Very robust and resistant



8305 // PROTECT 300

Material:	Nitrile, purple
Coating:	-
Size(s):	S-XL
Packaging unit:	10 piece(s)
Certification:	EN ISO 21420, EN ISO 374

NITRAS PROTECT 300, nitrile disposable gloves, purple (colour code: 6200), non-sterile, rolled edge, powder free, manufactured according to EN 455, for food contact, ambidextrous (left and right hand), medical gloves, micro-roughened fingertips, length approx. 300 mm, AQL 1.5, protection against micro-organisms, bacteria and viruses, particularly sensitive, perfect for working with small parts, very good tactile sensitivity, very good dry and oil grip, very robust and resistant, medical devices class I, PPE risk category III, box per 100 pieces

ISO 374-1/Type B



JKT

ISO 374-5:2016



VIRUS

EN ISO 21420

NI






AQL
1,5



TESTED
ACCORDING TO
EN 455

Industries

-  Cleaning Edition
-  Food Edition
-  Medical Edition

Advantages



- + Very good tactile sensitivity
- + Very good dry and oil grip
- + Very robust and resistant



Further images



Performance values

EN ISO 374-1:2016 + A1:2018	Test parameter	Performance levels	Test result
 JKT	J - n-Heptane	1-6	3 (> 60)
	K - Sodium hydroxide 40%	1-6	6 (> 480)
	T - Formaldehyde 37%	1-6	4 (> 120)
EN ISO 374-5:2016	Test parameter	Performance levels	Test result
 VIRUS	Penetration	-	Passed
	Protection against viruses	-	-
EN ISO 21420:2020	Test parameter	Performance levels	Test result
	Dexterity	1-5	5

Additional information

	Material	Colour
Palm	Nitrile	purple
Index finger and thumb	Nitrile	purple
Back of the hand	Nitrile	purple
Fingertips	Nitrile	purple
Cuff / knitted wrist	Nitrile	purple

Thickness

Finger	0,12 mm
Palm	0,08 mm
Cuff	0,08 mm

Size	Length (mm)	Width palm (mm)
S	300	84
M	300	84
L	300	105
XL	300	113

Properties	Before ageing	After ageing
Elongation at break	6%	6%

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
S	4059904096394	4059904096431	0,37 m	0,26 m	0,28 m	6,13 kg	70
M	4059904096400	4059904096448	0,37 m	0,26 m	0,28 m	6,26 kg	70
L	4059904096417	4059904096455	0,37 m	0,26 m	0,28 m	6,68 kg	70
XL	4059904096424	4059904096462	0,37 m	0,26 m	0,28 m	6,86 kg	70

Customs tariff number: 40151200000

Country of origin: Malaysia