

# 8321 // PINK WAVE

Nitrile disposable gloves  
PPE risk category III | Medical device class I



Very good tactile sensitivity  
Very good dry and oil grip  
Very robust and resistant



AQL  
1,5

NI

TESTED  
ACCORDING TO  
EN 455







## 8321 // PINK WAVE

<b>Material:</b>	Nitrile, light pink
<b>Coating:</b>	-
<b>Size(s):</b>	XS-L
<b>Packaging unit:</b>	10 piece(s)
<b>Certification:</b>	EN ISO 21420, EN ISO 374, EN 421

NITRAS PINK WAVE, nitrile disposable gloves, light pink (colour code: 6300), non-sterile, rolled edge, powder free, manufactured according to EN 455, ambidextrous (left and right hand), medical gloves, micro-roughened fingertips, colour code system for sizes, AQL 1.5, protection against micro-organisms, bacteria and viruses, protection against radioactive contamination by solid particles, particularly sensitive, perfect for working with small parts, very good tactile sensitivity, very good dry and oil grip, very robust and resistant, medical devices class I, PPE risk category III, box per 100 pieces



### Industries

-  Beauty Edition
-  Black Edition
-  Care Edition
-  Cleaning Edition
-  Dental Edition
-  Medical Edition

### Advantages

- + Very good tactile sensitivity
- + Very good dry and oil grip
- + Very robust and resistant



Further images

---



Colour variants

---



**8321 // PINK WAVE**  
Colour code: 6300



**8304 // GLAMOUR WAVE**  
Colour code: 6300



**8310 // WHITE WAVE**  
Colour code: 1100



**8311 // BLUE WAVE**  
Colour code: 2200



**8320 // BLACK WAVE**  
Colour code: 1000



**8322 // PURPLE WAVE**  
Colour code: 6200



**8323 // GREEN WAVE**  
Colour code: 3000



**8324 // INDIGO WAVE**  
Colour code: 2400



**8325 // MINT WAVE**  
Colour code: 3200



**8326 // RED WAVE**  
Colour code: 6000






**8327 // CANDY WAVE**  
Colour code: 6400



**8328 // VIOLET WAVE**  
Colour code: 6200

Performance values

EN ISO 374-1:2016 + A1:2018	Test parameter	Performance levels	Test result
 JKT	J - n-Heptane	1-6	3 (> 60)
	K - Sodium hydroxide 40%	1-6	6 (> 480)
	T - Formaldehyde 37%	1-6	4 (> 120)
EN ISO 374-5:2016	Test parameter	Performance levels	Test result
 VIRUS	Penetration	-	Passed
	Protection against viruses	-	Passed
EN ISO 21420:2020 + A1:2024	Test parameter	Performance levels	Test result
	Dexterity	1-5	5
EN 421:2010	Test parameter	Performance levels	Test result
			

**Additional information**

	Material	Colour
Palm	Nitrile	light pink
Index finger and thumb	Nitrile	light pink
Back of the hand	Nitrile	light pink
Fingertips	Nitrile	light pink
Cuff / knitted wrist	Nitrile	light pink

**Thickness**

Finger	0,10 mm
Palm	0,08 mm
Cuff	0,07 mm

Size	Length (mm)	Width palm (mm)
XS	240	76
S	240	84
M	240	94
L	240	
XL	240	113

Properties	Before ageing	After ageing
Elongation at break	500%	400%
Tear resistance	14 MPa	14 MPa

**Logistics data**

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
XS	4059904035584	4059904003996	0,37 m	0,23 m	0,22 m	3,95 kg	80
S	4059904035591	4059904004009	0,37 m	0,23 m	0,22 m	4,38 kg	80
M	4059904035607	4059904004016	0,37 m	0,23 m	0,22 m	4,56 kg	80
L	4059904035614	4059904004023	0,37 m	0,23 m	0,22 m	5,07 kg	80
XL	4059904035621	4059904004030	0,37 m	0,23 m	0,22 m	5,58 kg	80

Customs tariff number: 40151200000

Country of origin: Malaysia