

8322 // PURPLE WAVE

Nitrile disposable gloves
PPE risk category III | Medical device class I




Very good tactile sensitivity
Very good dry and oil grip
Very robust and resistant



8322 // PURPLE WAVE

Material:	Nitrile, purple
Coating:	-
Size(s):	XS-XL
Packaging unit:	10 piece(s)
Certification:	EN ISO 21420, EN ISO 374, EN 421

NITRAS PURPLE WAVE, nitrile disposable gloves, purple (colour code: 6200), non-sterile, rolled edge, powder free, manufactured according to EN 455, for food contact, ambidextrous (left and right hand), medical gloves, micro-roughened fingertips, colour code system for sizes, AQL 1.5, protection against micro-organisms, bacteria and viruses, protection against radioactive contamination by solid particles, particularly sensitive, perfect for working with small parts, very good tactile sensitivity, very good dry and oil grip, very robust and resistant, medical devices class I, PPE risk category III, box per 100 pieces



Industries

-  Beauty Edition
-  Black Edition
-  Cleaning Edition
-  Dental Edition
-  Food Edition
-  Medical Edition

Advantages

- + Very good tactile sensitivity
- + Very good dry and oil grip
- + Very robust and resistant



Further images



Colour variants



8322 // PURPLE WAVE
Colour code: 6200



8304 // GLAMOUR WAVE
Colour code: 6300



8310 // WHITE WAVE
Colour code: 1100



8311 // BLUE WAVE
Colour code: 2200



8320 // BLACK WAVE
Colour code: 1000



8321 // PINK WAVE
Colour code: 6300



8323 // GREEN WAVE
Colour code: 3000



8324 // INDIGO WAVE
Colour code: 2400



8325 // MINT WAVE
Colour code: 3200



8326 // RED WAVE
Colour code: 6000



8327 // CANDY WAVE
Colour code: 6400



8328 // VIOLET WAVE
Colour code: 6200

Performance values

EN ISO 374-1:2016 + A1:2018	Test parameter	Performance levels	Test result
 JKT	J - n-Heptane	1-6	3 (> 60)
	K - Sodium hydroxide 40%	1-6	6 (> 480)
	T - Formaldehyde 37%	1-6	4 (> 120)
EN ISO 374-5:2016	Test parameter	Performance levels	Test result
 VIRUS	Penetration	-	Passed
	Protection against viruses	-	Passed
EN ISO 21420:2020 + A1:2024	Test parameter	Performance levels	Test result
	Dexterity	1-5	5
EN 421:2010	Test parameter	Performance levels	Test result
			

Additional information

	Material	Colour
Palm	Nitrile	purple
Index finger and thumb	Nitrile	purple
Back of the hand	Nitrile	purple
Fingertips	Nitrile	purple
Cuff / knitted wrist	Nitrile	purple

Thickness (approx.)

Finger	0,10 mm
Palm	0,08 mm
Cuff	0,07 mm

Size	Length (mm)	Width palm (mm)
XS	240	76
S	240	84
M	240	94
L	240	105
XL	240	113

Properties	Before ageing	After ageing
Elongation at break	500%	400%
Tear resistance	14 MPa	14 MPa

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
XS	4059904035638	4059904004047	0,37 m	0,24 m	0,22 m	3,56 kg	80
S	4059904035645	4059904004054	0,37 m	0,24 m	0,22 m	4,03 kg	80
M	4059904035652	4059904004061	0,37 m	0,24 m	0,22 m	4,35 kg	80
L	4059904035669	4059904004078	0,37 m	0,24 m	0,22 m	4,44 kg	80
XL	4059904035676	4059904004085	0,37 m	0,24 m	0,22 m	4,85 kg	80

Customs tariff number: 40151200000

Country of origin: Malaysia