

8323 // GREEN WAVE

Nitrile disposable gloves
PPE risk category III | Medical device class I



Very durable
Micro-roughened finger tips
Very good dry and oil grip



8323 // GREEN WAVE

Material:	Nitrile, green
Coating:	-
Size(s):	XS-XL
Packaging unit:	10 piece(s)
Certification:	EN ISO 21420, EN ISO 374, EN 421

NITRAS GREEN WAVE, Nitrile disposable gloves, green (color code: 3000), non-sterile, rolled cuff, powder-free, manufactured according to EN 455, for food contact, ambidextrous (left and right hand), medical gloves, micro-textured fingertips, color-coded size system, AQL 1.5, protection against microorganisms, bacteria, and viruses, protection against radioactive contamination by solid particles, highly sensitive, ideal for working with small parts, excellent tactile sensitivity, very good dry and oily grip, particularly robust and durable, Class I medical devices, PPE risk category III, box of 100

ISO 374-1/Type B



JKT

ISO 374-5:2016



VIRUS

EN ISO 21420

EN 421:2010



NI

AQL
1,5

TESTED
ACCORDING TO
EN 455

Industries

-  Beauty Edition
-  Black Edition
-  Cleaning Edition
-  Dental Edition
-  Food Edition
-  Medical Edition

Advantages

- + Very durable
- + Micro-roughened finger tips
- + Very good dry and oil grip



Further images



Colour variants



8323 // GREEN WAVE
Colour code: 3000



8304 // GLAMOUR WAVE
Colour code: 6300



8310 // WHITE WAVE
Colour code: 1100



8311 // BLUE WAVE
Colour code: 2200



8320 // BLACK WAVE
Colour code: 1000



8321 // PINK WAVE
Colour code: 6300



8322 // PURPLE WAVE
Colour code: 6200



8324 // INDIGO WAVE
Colour code: 2400



8325 // MINT WAVE
Colour code: 3200



8326 // RED WAVE
Colour code: 6000



8327 // CANDY WAVE
Colour code: 6400



8328 // VIOLET WAVE
Colour code: 6200

Performance values

EN ISO 374-1:2016 + A1:2018	Test parameter	Performance levels	Test result
 JKT	J - n-Heptane	1-6	2 (> 30)
	K - Sodium hydroxide 40%	1-6	6 (> 480)
	T - Formaldehyde 37%	1-6	4 (> 120)
EN ISO 374-5:2016	Test parameter	Performance levels	Test result
 VIRUS	Penetration	-	Passed
	Protection against viruses	-	Passed
EN ISO 21420:2020 + A1:2024	Test parameter	Performance levels	Test result
	Dexterity	1-5	5
EN 421:2010	Test parameter	Performance levels	Test result
			

Additional information

	Material	Colour
Palm	Nitrile	green
Index finger and thumb	Nitrile	green
Back of the hand	Nitrile	green
Fingertips	Nitrile	green
Cuff / knitted wrist	Nitrile	green

Thickness (approx.)

Finger	0,10 mm
Palm	0,08 mm
Cuff	0,07 mm

Size	Length (mm)	Width palm (mm)
XS	240	76
S	240	84
M	240	94
L	240	105
XL	240	113

Properties	Before ageing	After ageing
Elongation at break	500%	400%
Tear resistance	6 MPa	14 MPa

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
XS	4059904230712	4059904230705	0,37 m	0,24 m	0,22 m	3,56 kg	80
S	4059904230736	4059904230729	0,37 m	0,24 m	0,22 m	4,03 kg	80
M	4059904230750	4059904230743	0,37 m	0,24 m	0,22 m	4,35 kg	80
L	4059904230774	4059904230767	0,37 m	0,24 m	0,22 m	4,4 kg	80
XL	4059904230798	4059904230781	0,37 m	0,24 m	0,22 m	4,85 kg	80

Customs tariff number: 40151200000

Country of origin: Malaysia