

8324 // INDIGO WAVE

Nitrile disposable gloves
PPE risk category III | Medical device class I



Very good tactile sensitivity
Very good dry and oil grip
Very robust and resistant



AQL
1,5

NI

TESTED
ACCORDING TO
EN 455

8324 // INDIGO WAVE

Material:	Nitrile, dark blue
Coating:	-
Size(s):	XS-XL
Packaging unit:	10 piece(s)
Certification:	EN ISO 21420, EN ISO 374, EN 421

NITRAS INDIGO WAVE, nitrile disposable gloves, dark blue (colour code: 2400), non-sterile, rolled edge, powder free, manufactured according to EN 455, ambidextrous (left and right hand), medical gloves, micro-roughened fingertips, colour code system for sizes, AQL 1.5, protection against micro-organisms, bacteria and viruses, protection against radioactive contamination by solid particles, particularly sensitive, perfect for working with small parts, very good tactile sensitivity, very good dry and oil grip, very robust and resistant, medical devices class I, PPE risk category III, box per 100 pieces

ISO 374-1/Type B



JKT

ISO 374-5:2016



VIRUS

EN ISO 21420

EN 421:2010


AQL
1,5

NI

TESTED
ACCORDING TO
EN 455

Industries

-  Beauty Edition
-  Black Edition
-  Care Edition
-  Cleaning Edition
-  Dental Edition
-  Medical Edition

Advantages

- + Very good tactile sensitivity
- + Very good dry and oil grip
- + Very robust and resistant



Further images



Colour variants



8324 // INDIGO WAVE
Colour code: 2400



8304 // GLAMOUR WAVE
Colour code: 6300



8310 // WHITE WAVE
Colour code: 1100



8311 // BLUE WAVE
Colour code: 2200



8320 // BLACK WAVE
Colour code: 1000



8321 // PINK WAVE
Colour code: 6300



8322 // PURPLE WAVE
Colour code: 6200



8323 // GREEN WAVE
Colour code: 3000



8325 // MINT WAVE
Colour code: 3200



8326 // RED WAVE
Colour code: 6000



8327 // CANDY WAVE
Colour code: 6400



8328 // VIOLET WAVE
Colour code: 6200

Performance values

EN ISO 374-1:2016 + A1:2018	Test parameter	Performance levels	Test result
 JKT	J - n-Heptane	1-6	3 (> 60)
	K - Sodium hydroxide 40%	1-6	6 (> 480)
	T - Formaldehyde 37%	1-6	4 (> 120)
EN ISO 374-5:2016	Test parameter	Performance levels	Test result
 VIRUS	Penetration	-	Passed
	Protection against viruses	-	Passed
EN ISO 21420:2020 + A1:2024	Test parameter	Performance levels	Test result
	Dexterity	1-5	5
EN 421:2010	Test parameter	Performance levels	Test result
			

Additional information

	Material	Colour
Palm	Nitrile	dark blue
Index finger and thumb	Nitrile	dark blue
Back of the hand	Nitrile	dark blue
Fingertips	Nitrile	dark blue
Cuff / knitted wrist	Nitrile	dark blue

Thickness (approx.)

Finger	0,10 mm
Palm	0,08 mm
Cuff	0,07 mm

Size	Length (mm)	Width palm (mm)
XS	240	76
S	240	84
M	240	94
L	240	105
XL	240	113

Properties	Before ageing	After ageing
Elongation at break	500%	400%
Tear resistance	14 MPa	14 MPa

Logistics data

Size	GTIN (SPU)	GTIN (PU)	Length (m)	Width (m)	Height (m)	Weight (kg)	Quantity PU / pallet
XS	4059904036826	4059904005235	0,37 m	0,24 m	0,22 m	3,49 kg	80
S	4059904036833	4059904005242	0,37 m	0,24 m	0,22 m	3,97 kg	80
M	4059904036840	4059904005259	0,37 m	0,24 m	0,22 m	4,45 kg	80
L	4059904036857	4059904005266	0,37 m	0,24 m	0,22 m	4,16 kg	80
XL	4059904036864	4059904005273	0,37 m	0,24 m	0,22 m	4,8 kg	80

Customs tariff number: 40151200000

Country of origin: Malaysia